

Covent Test Plan

SR01	Mandatory Ventilation							
SR08	PEEP Pressure 5, 10, 15 cmH2O in increments of 5 cmH2O (+/- 5 cmH2O) (see volume control testing)							
SR09	Mandatory Ventilation - I:E 1:2, 1:3, and 1:4 options available **1:5 & 1:1 I:E per respiratory rate requirements (see I:E column of Volume Control testing)							
SR10	Respiratory Rate 10-30 breaths per minute in increments of 2 bpm 10, 12, 15, 20, 30 to follow I:E ratios (see volume control testing or demonstrate individually ?)							
SR11	Tidal Volume Option #1: Input height and gender for 6cc/kg TV (+/-10% or 10mL) Option #2: 350cc (for average woman) and 450cc (for average man) (+/-10% or 10mL) Option #3: 400cc only (+/-10% or 10mL) Option #4: 300, 400, 500, 600 (+/-10% or 100mL) **200 ISO only (see TidalVolume column in Volume Control Testing)							
SR02	Spontaneous Ventilation see spontaneous.mov							
SR06	Apnea back up kicks in at 30 or 60 seconds (+/- 5sec) see spontaneous.mov							
SR03	Volume Control							
Section 201.12.1.101 from ISO 80601-2-12								
	Compliance ml/hPa ±10%	Resistance hPa/l/s ±10%	Tidal Volume ml	Breaths / min (a)	Inspire time (s)	O2	BAP (PEEP)	I:E
1	50	5	500	20	1	30	5	1:2
	Pressure Difference : -1.14% Tidal Volume Difference : +32.60%							

2	50	20	500	12	1	90	10	1:4
	Pressure Difference : +4.31% Tidal Volume Difference : +57.49%							
3	20	5	500	20	1	90	5	1:2
	Pressure Difference : -0.27% Tidal Volume Difference : +56.02%							
4	20	20	500	20	1	30	10	
5	20	20	300	20	1	30	5	
6	20	50	300	12	1	90	10	1:3
7	10	50	300	20	1	30	10	
8	10	10	200	20	1	90	5	
C200	50	5	200	15	1	20	5	1:3
	Pressure Difference : +2.15% Tidal Volume Difference : +67.67%							
C300	50	5	300	15	1	20	5	1:3
C400	50	5	400	15	1	20	5	1:3
C500	50	5	500	15	1	20	5	1:3
C600	50	5	600	15	1	20	5	1:3
P15	50	5	500	10	1	20	15	1:5
	Pressure Difference : -14.44% Tidal Volume Difference : -11.64%							
R30	50	5	200	30	1	20	5	1:1
	Pressure Difference : +2.55% Tidal Volume Difference : +46.53%							
a) if the end expiratory flow does not reach zero, reduce set rate until it does								
SR04	Pressure Control 5-60 +/- 5 cmH2O Control Module							
	Compliance ml/hPa ±10%	Resistance hPa//s ±10%	Tidal Volume ml (a)	Breaths / min (b)	Inspire time (s) (c)	O2	BAP (PEEP)	Δ insp pressure (d)
1	50	5	500	20	1	30	5	10

	Pressure Difference : -2.37% Tidal Volume Difference : +17.60%							
2	50	20	500	12	1	90	10	15
	Pressure Difference : -0.09% Tidal Volume Difference : +31.01%							
3	20	5	500	20	1	90	5	25
	a) <i>intended tidal volume for the selection of test conditions on test lung</i> b) <i>if the end expiratory flow does not reach zero, reduce set rate until it does</i> c) <i>the rise time of the ventilator should be set to a value that insecure pressure can be reached</i> d) <i>for this test the set pressure is relative to BAP (PEEP)</i>							
SR05	Pressure Support 10-15 +/- 5 cmH2O flow or pressure triggered							
SR07	Flow Rate > 60 liters per minute							
C600	refer to C600 test data							
Alarms								
SR19	air_alarm							
SR20	pressure_alarm							
SR21	spontaneous							
SR22	peep_alarm							
SR23	air_alarm							
SR24	air_alarm							
SR25	decibel.MP4							